

## Waive™ Test Cup

- Two drugs per test strip
- Operator activated – activate test when ready
- Detects twice as many BZO's as other CLIA-waved cups
- Leak-proof urine collection and testing system
- No maximum specimen volume
- No donor access to test
- Neat, clean, extremely easy
- No urine exposure
- Detects up to 12 drugs and adulterants



**STEP 1** *Secure test lid to specimen cup*



**STEP 2** *Lay cup down to activate test*



**STEP 3** *Interpret results*

**NOTE:** Do not tip the cup upside down. To activate, lay the cup on its side on the tabs.