

Now test for all the 'usual suspects' and alcohol

with

DRUGCHECK[®]
NxStep OnSite

Shown with adulterant and alcohol test strips



NEW!
1st Urine Alcohol Test Cup!

Available with test strips for alcohol and adulterants, or alcohol only



Alcohol only

NEW!

1st Urine Alcohol Test Cup!

DRUGCHECK[®]

NxStep OnSite

15 tests in one



- ✓ 15-panel test
- ✓ Available with alcohol and/or adulterant tests
- ✓ 8 oz. cup
- ✓ 9% larger opening
- ✓ Unbreakable cup surface

- ✓ Leak-proof transport technology
- ✓ 15 individual chambers eliminate drug interpretation confusion
- ✓ Tamper-proof strip insert
- ✓ Made in U.S.A.
- ✓ Results in 5 minutes



More great innovations from the leader


The manufacturer of the DrugCheck[®] Drug Test Cup and Dip Drug Test proudly introduces the first 15-panel drug test cup on the market. With the new NxStep OnSite Drug Test Cup, you can now test for up to 15 drug types in a single cup with **results in 5 minutes**.

- ✓ NO ASSEMBLY
- ✓ NO TILTING TO ACTIVATE
- ✓ NO SPLIT SAMPLES OR CHAMBER TO RELY ON
- ✓ NO CARD DIPPING
- ✓ NO CROSS CONTAMINATION
- ✓ NO MOVING PARTS

Designed for the person giving the test

With its superior design, the DrugCheck[®] NxStep OnSite Drug Test Cup requires no handling or manipulation to activate testing. The NxStep OnSite is **FDA cleared** and has the same accuracy when compared to certified lab testing.

AMP AMP AMP AMP	Amphetamine	1000 ng/ml
BAR BAR BAR BAR	Barbiturate	300 ng/ml
BZD BZD BZD	Benzodiazepine	300 ng/ml
COC COC COC COC	Cocaine	300 ng/ml
MET MET MET MET	Methamphetamine	1000 ng/ml
MTD MTD MTD	Methadone	300 ng/ml
OPI OPI OPI OPI	Opiates	300/2000 ng/ml
OXY OXY OXY OXY	Oxycodone	100 ng/ml
PCP PCP PCP PCP	Phencyclidine (PCP)	25 ng/ml
THC THC THC	Marijuana	50 ng/ml
MDMA MDMA MDMA	Ecstasy	500 ng/ml
PPX PPX PPX	Propoxyphene	300 ng/ml
TCA TCA TCA	Tricyclic Antidepressants	1000 ng/ml



For use with alcohol test strips

For use with adulterant test strips

	ABNORMAL (LOW)	NORMAL	ABNORMAL (HIGH)
pH Immediate	2 3	4 5 6 7 8	>10
Specific Gravity (SG) 45 seconds	1.00	1.005 1.015 1.025	>1.030
Nitrite (NI) 45 seconds	0	0.1-0.2 0.3-0.4	>0.4
Creatinine (CR) 45 seconds	Negative	10 20 30 100 200 mg/dl	

Distributed by:
Preferred Drug Testing
 PO Box 14 - Worthington, MN. 56187
 605-366-2994
 email: info@drugcup.com
www.drugcup.com