



Collect sufficient sample the first time with the first ever sponge saturation indicator!

- Detect up to 10 drugs
- 100% tamper evident
- Saliva held in airtight chamber for confirmation
- Photocopy or scan results
- Made in the USA

Oral Fluid Collection Media (sponge)

- Highest recovery of DOA analytes
- Buffer not required
- Highest absorption rate
- Reduces viscosity

SPONGE SATURATION INDICATOR



FORENSIC USE ONLY

SALIVA CUTOFFS	
Alcohol (ALC)	.02%
Amphetamine (AMP)	50
Benzodiazepines (BZO)	10
Buprenorphine (BUP)	5
Cocaine (COC)	20
Cotinine (COT)	50
EDDP	20
Ketamine (KET)	50
Marijuana (THC)	12 (parent) 50
Methadone (MTD)	30
Methamphetamine (MET)	50
Opiates (OPI)	40
Oxycodone (OXY)	40
Phencyclidine (PCP)	10

All cutoffs in ng/mL except alcohol

Distributed by:

Preferred Drug Testing • PO Box 14 Worthington, Minnesota 56187 • 605-366-2994 • www.drugcup.com